MiniCLEAN Update

$\bullet \bullet \bullet$

Nathan Eggen

Previously at MiniCLEAN



- Successfully pumped and purged detector while warming up cold fingers.
- Had confusion with Argon billing to PNNL
 - Had none to fill with over 6 days with no underground access.
- Current IV temp is ~133 K at the bottom and ~135 at the top

Updating rat:

- Overall goal to upgrade MiniCLEAN version of RAT to use Geant 4.10
- Met with Hans to help him get started with RAT and using HEP cluster.